

Veriport API Documentation

Automated Report Delivery via TPA API (SNS Integration)

Overview

We have implemented a system that automatically delivers lab reports to a Third-Party Administrator (TPA) through Amazon SNS (Simple Notification Service).

- Reports are serialized into XML files
- Securely stored in Amazon Web Services (AWS) Simple Storage Service (S3)
- Made available for download via a one-time use, unique, short URL
- The TPA receives a notification containing URL via AWS Simple Notification Service (SNS)
<https://aws.amazon.com/sns/>

👉 This ensures secure, trackable, and automated report delivery.

Workflow

When a report is generated, the system performs the following steps:

1. Generate XML file
 - The report is serialized into XML format.
2. Upload to S3
 - The XML file is uploaded to a secure S3 bucket.
3. Generate presigned URL
 - A presigned URL is created for the uploaded file.
 - Allows temporary, secure access to the file without exposing S3 credentials.
4. Wrap into a short URL
 - The presigned URL is wrapped in a short URL that points to our server.
 - This allows us to register access before redirecting to the actual S3 file.
5. Publish SNS notification
 - A message is sent to a dedicated SNS topic associated with the TPA.
 - The message contains the short URL, which the TPA can use to download the XML report.

TPA Requirements

To receive reports, the TPA must provide the following:

1. HTTPS Endpoint for SNS messages

The TPA must expose an HTTPS endpoint that can handle SNS messages. Example URL:
<https://api.tpa.com/sns-receiver>

Sample Django REST Framework Implementation:

```
import json
import requests
from rest_framework.views import APIView
from rest_framework.response import Response
from rest_framework import permissions

class SNSReceiverView(APIView):
    """
    Example implementation for receiving SNS notifications.
    """
    permission_classes = [permissions.AllowAny]

    def post(self, request, format=None):
        try:
            body = json.loads(request.body)
            message_type = body.get('Type')
            message = body.get('Message')

            if message_type == 'SubscriptionConfirmation':
                subscribe_url = body.get('SubscribeURL')
                if subscribe_url:
                    response = requests.get(subscribe_url)
                    if response.status_code == 200:
                        return Response({'message': 'Subscription confirmed'}, status=200)
                    return Response({'message': 'Failed to confirm subscription'}, status=400)

            if message_type == 'Notification':
                notification_message = json.loads(message)
                download_url = notification_message.get('download_url')
                print(f"Download file here: {download_url}")
                # Implement file download & processing here

            return Response({'message': 'Message processed'}, status=200)
        except json.JSONDecodeError:
            return Response({'message': 'Invalid JSON'}, status=400)
```

 The endpoint must be publicly accessible via HTTPS so AWS SNS can post notifications.

2. Enable Auto-Delivery Setting in Veriport

On the TPA side, the auto-delivery account setting must be enabled to allow auto-delivery of MRO reports to an external system.

Manual Configuration (Performed by Us)

Once the TPA provides their endpoint:

1. Create an SNS topic in AWS for the TPA account.
 2. Add the provided HTTPS receiver URL as a subscription.
 3. Wait for the TPA to confirm the subscription (via [SubscriptionConfirmation](#)).
 4. Once confirmed, and the auto-delivery flag is enabled on their account, the TPA will begin receiving messages containing report download URLs.
-

Summary

- Reports are serialized to XML, uploaded to S3, and delivered via SNS.
 - TPAs receive notifications with short URLs for secure, trackable downloads.
 - The TPA must:
 - Provide a valid HTTPS receiver endpoint
 - Enable auto-delivery in account settings
 - We handle SNS topic creation, subscription, and delivery configuration on our side.
-

Schema 1.5:

Variable Name	Expected Value	Explanation
VeriportReportId	56789987	Unique Veriport Id, updates to data will be based on this Id (numeric).
ReportRevisionNumber	rev1	New versions of the MRO report will have incremental revisions (rev1, rev2, etc.). The first report may have null.
ReportUpdate	false	Indicates if changes are made that are not material to the actual MRO report (e.g., transmitting the SignedCcf2).
ReasonForTestLab	RAND	Reason for test code from lab.
ReasonForTest	RAN	Veriport database value for reason for test based on TPA account settings.

ReasonForTestName	Random	Reason for test mapped description (Return To Duty, Reasonable Suspicion / Cause, Random, etc.).
Regulated	true	Indicates whether the test is regulated (e.g., by DOT).
TestingAuthority	DOT	Governing authority overseeing the test; HHS, NRC, DOT.
DotAgency	FMCSA	Specific DOT agency; FAA, FMCSA, FRA, FTA, PHMSA, USCG.
CustomerAccountId	1234	Unique Veriport Account ID for the customer/parent (TPA) account (numeric).
CustomerName	ABC Testing Services Inc.	Name of the customer/parent account.
CustomerAddress1	123 Main St.	First line of the customer's address.

CustomerAddress2	Unit D	Second line of the customer's address.
CustomerCity	Barrie	City where the customer is located.
CustomerStateProv	ON	Province/state (2-character standard).
CustomerZipPostalCode	L4N 8W8	Postal code of the customer.
CustomerCountry	CA	Country where the customer is located; CA, US, MX.
CustomerPhone	705-444-2121	Customer phone number in xxx-xxx-xxxxx format.
ClientAccountId	4321	Unique Veriport Account ID for the client/child (Employer) account (numeric).
ClientCustomerAccountNumber	B1234	Customer-assigned identifier for client/child (Employer) account (alphanumeric).

ClientName	XYZ Trucking Inc.	Employer name.
ClientAddress1	123 Main St.	First line of employer address.
ClientAddress2	Unit D	Second line of employer address.
ClientCity	Barrie	City where employer is located.
ClientStateProv	ON	Province/state (2-character standard).
ClientZipPostalCode	L4N 8W8	Postal code of employer.
ClientCountry	CA	Employer country.
ClientPhone	705-444-2121	Employer phone number.
ClientNonRegQuant	true	Indicates whether non-regulated quantitative results should be included.
DonorSystemId	D4R7TY	Unique Veriport Donor System Id (alphanumeric).

DonorFirst	Michael	Donor's first name.
DonorMiddle	Paul	Donor's middle name.
DonorLast	Smith	Donor's last name.
DonorPhone	705-333-1212	Donor phone number.
DonorEmail	mps@email.com	Donor email.
DonorId	H78975681139	Donor's identifier (e.g., driver's license).
DonorIdType	DL	Type of Id provided by donor (DL, Passport, etc.).
DonorIdCdIStateProv	ON	Province/state for donor's DL.
DonorIdCdICountry	CA	Country for donor's DL.
DonorDob	05-23-1990	Donor's date of birth (MM-DD-YYYY).
MroAccountId	8798	Unique Veriport Account Id for the MRO.

MroCompanyName	Super MRO	Name of the MRO company.
MroAddress1	123 Main St	First line of MRO's address.
MroAddress2	Suite 1000	Second line of MRO's address.
MroCity	Houston	City where MRO is located.
MroStateProv	TX	Province/state of MRO.
MroZipPostalCode	77057	Postal code of MRO.
MroCountry	US	Country of MRO.
MroPhone	705-333-1212	MRO phone number.
MroFax	705-222-1212	MRO fax number.
MroEmail	support@supermro.com	MRO email.
CollectionSiteId	34567	Unique Veriport Site Id (numeric).

CollectionSiteCode	Site 1234	Collection site code from lab (alphanumeric).
CollectionSiteCustomerCode	Site 4321	Collection site code from customer (alphanumeric).
CollectionSiteName	Main Office	Name of collection site.
CollectionSiteAddress1	123 Main St.	First line of collection site address.
CollectionSiteAddress2	Unit D	Second line of collection site address.
CollectionSiteCity	Barrie	Collection site city.
CollectionSiteStateProv	ON	Province/state of collection site.
CollectionSiteZipPostalCode	L4N 8W8	Postal code of collection site.
CollectionSiteCountry	CA	Country of collection site.
CollectionSitePhone	705-444-2121	Collection site phone.

CollectorName	Jennifer Smith	Collector's name.
CollectorPhone	705-555-1212	Collector's phone number.
CollectorFax	705-555-2121	Collector's fax number.
CollectorOther	alpha/numeric	Text field as entered on CCF.
CollectorEmail	collector@collection.com	Collector's email (from lab/eCCF).
LabName	Clinical Reference Laboratory, Inc.	Name of laboratory.
LabId	45	Unique Id for lab in Veriport system.
LabAddress1	123 Main St.	First line of lab address.
LabAddress2	Unit D	Second line of lab address.
LabCity	Barrie	City of lab.
LabStateProv	ON	Province/state of lab.

LabZipPostalCode	L4N 8W8	Postal code of lab.
LabCountry	CA	Country of lab.
LabPhone	705-444-2121	Lab phone.
LabFax	705-444-1212	Lab fax.
LabAccountNumber	45678	Lab account number.
LabSubaccountNumber	9876	Lab sub-account number.
LabSpecimenId	CF18142152	Specimen Id (not unique, do not rely as reference).
LabAccessionNumber	12345678	Lab accession number.
LabCustomerReferenceId	AU123456	Reference number provided to lab by customer (per test basis).
LabCertifyingScientist	Brianna Dingman	Scientist certifying lab results.
MroName	Dr. Alfred Smith	Name of Medical Review Officer.

MroQualifyingOrganization	MROCC	Qualifying organization for MRO.
MroRegistrationNumber	1234567890	MRO registration number with qualifying organization.
MroOverallResult	Positive Dilute	MRO-defined overall result.
MroMisOverallResult	Positive	Mapped DOT MIS result value (Positive, Negative, Refusal, Cancelled, etc.).
LabPanelCode	W215	Code for the lab panel performed.
LabPanelName	Urine Drug Test	Name of the lab test panel.
LabSpecimenType	Urine	Specimen type tested.
LabTestName / TestCodeName	e.g., Cocaine	Name of specific drug test.
LabTestCode / TestCodeNumber	e.g., COC	Code of the drug test.
UnitOfMeasure	ng/mL	Measurement unit.

ScreenCutoff	varies by drug	Screening cutoff value.
ConfirmCutoff	varies by drug	Confirmation cutoff value.
MroResult	Positive/Negative	Result for specific drug/panel.
NonRegQuant	(if applicable)	Numeric value if non-regulated & allowed.
SpecimenCollectedDate	05-23-2024	Date specimen collected (MM-DD-YYYY).
MroReceivedCcf2Date	05-23-2024	Date MRO received CCF.
LabReceivedDate	05-24-2024	Date lab received specimen.
LabReceivedTime	17:05:00Z	Time lab received specimen (HH:mm:ssZ).
LabResultedToMroDate	05-26-2024	Date lab reported results to MRO.

LabResultedToMroTime	10:09:05Z	Time lab reported results to MRO (HH:mm:ssZ).
MroReportedDate	05-26-2024	Date MRO reported final results.
MroReportedTime	16:30:09Z	Time MRO reported final results.
MroRecalledDate	(nullable)	Date when MRO recalled report, if applicable.
MroRecalledTime	(nullable)	Time when MRO recalled report.
MroRereleasedDate	(nullable)	Date MRO re-released report after recall.
MroRereleasedTime	(nullable)	Time MRO re-released report.
MroLetterComment	text	MRO letter comment.
MroSafetyAdvisoryComment	text	MRO safety advisory comment.
MroLetter	Base64 Encoded PDF	MRO letter.

SignedCcf2	Base64 Encoded JPG	The copy 2 of the CCF will be signed at a later date and then an updated file will be available to import with this new document.
Cg719P	Base64 Encoded JPG	Coast Guard 719P form, when applicable

Sample XML:

```
<?xml version="1.0" encoding="UTF-8"?>
<VeriportMroData version="1.5">
  <VeriportReportId>56789987</VeriportReportId> <!-- IMPORTANT Unique Veriport Id, updates
to data will be based on the VeriportReportId - numeric -->
  <TestInformation>
    <ReportRevisionNumber>rev1</ReportRevisionNumber> <!-- new versions of the MRO report
will have incremental revisions, rev1, rev2, etc. The first report will have a null value
-->
    <ReportUpdate>false</ReportUpdate> <!-- if changes are made to the report information,
that are not material to the actual MRO report (therefore not revised), this value will be
"true" e.g. transmitting the SignedCcf2 -->
    <ReasonForTestLab>RAND</ReasonForTestLab> <!-- Reason for test code from lab -->
    <ReasonForTest>RAN</ReasonForTest> <!-- Veriport database value for reason for test
(does not change) -->
    <ReasonForTestName>Random</ReasonForTestName> <!-- Reason for test mapped description
Return To Duty, Reasonable Suspicion / Cause, Random, Pre-Employment, Post-Accident,
Follow-up, Other (may change)-->
    <Regulated>true</Regulated>
    <TestingAuthority>DOT</TestingAuthority> <!-- HHS NRC DOT -->
    <DotAgency>FMCSA</DotAgency> <!-- FAA FMCSA FRA FTA PHMSA USCG -->
  </TestInformation>
  <Customer>
    <CustomerAccountId>1234</CustomerAccountId> <!-- Unique Veriport Account ID for the
customer/parent (TPA) account - numeric -->
    <CustomerName>ABC Testing Services Inc.</CustomerName>
    <CustomerAddress1>123 Main St.</CustomerAddress1>
    <CustomerAddress2>Unit D</CustomerAddress2>
    <CustomerCity>Barrie</CustomerCity>
    <CustomerStateProv>ON</CustomerStateProv>
    <CustomerZipPostalCode>L4N 8W8</CustomerZipPostalCode>
    <CustomerCountry>CA</CustomerCountry>
    <CustomerPhone>705-444-2121</CustomerPhone> <!-- phone in xxx-xxx-xxxxx format always
-->
  </Customer>
  <Client>
    <ClientAccountId>4321</ClientAccountId> <!-- Unique Veriport Account ID for the
client/child (Employer) account - numeric -->
    <ClientCustomerAccountNumber>B1234</ClientCustomerAccountNumber> <!-- Customer assigned
for the client/child (Employer) account - alpha/numeric -->
    <ClientName>XYZ Trucking Inc.</ClientName>
    <ClientAddress1>123 Main St.</ClientAddress1>
    <ClientAddress2>Unit D</ClientAddress2>
    <ClientCity>Barrie</ClientCity>
    <ClientStateProv>ON</ClientStateProv>
    <ClientZipPostalCode>L4N 8W8</ClientZipPostalCode>
    <ClientCountry>CA</ClientCountry>
    <ClientPhone>705-444-2121</ClientPhone>
    <ClientNonRegQuant>true</ClientNonRegQuant>
```

```

</Client>
<Donor>
  <DonorSystemId>D4R7TY</DonorSystemId> <!-- Unique Veriport Donor System Id -
alphanumeric-->
  <DonorFirst>Michael</DonorFirst>
  <DonorMiddle>Paul</DonorMiddle>
  <DonorLast>Smith</DonorLast>
  <DonorPhone>705-333-1212</DonorPhone>
  <DonorEmail>mps@email.com</DonorEmail>
  <DonorId>H78975681139</DonorId>
  <DonorIdType>DL</DonorIdType>
  <DonorIdCdlStateProv>ON</DonorIdCdlStateProv> <!-- if "DonorIdType" is "DL" -->
  <DonorIdCdlCountry>CA</DonorIdCdlCountry> <!-- if "DonorIdType" is "DL" -->
  <DonorDob>05-23-1990</DonorDob> <!-- Dates in mm-dd-yyyy format -->
</Donor>
<Mro>
  <MroAccountId>8798</MroAccountId> <!-- Unique Veriport Account Id -->
  <MroCompanyName>Super MRO</MroCompanyName>
  <MroAddress1>123 Main St</MroAddress1>
  <MroAddress2>Suite 1000</MroAddress2>
  <MroCity>Houston</MroCity>
  <MroStateProv>TX</MroStateProv>
  <MroZipPostalCode>77057</MroZipPostalCode>
  <MroCountry>US</MroCountry>
  <MroPhone>705-333-1212</MroPhone>
  <MroFax>705-222-1212</MroFax>
  <MroEmail>support@supermro.com</MroEmail>
</Mro>
<CollectionSite>
  <CollectionSiteId>34567</CollectionSiteId> <!-- Unique Veriport Site Id - numeric-->
  <CollectionSiteCode>Site 1234</CollectionSiteCode> <!-- Collection Site code as provided
by the lab - alpha/numeric -->
  <CollectionSiteCustomerCode>Site 4321</CollectionSiteCustomerCode> <!-- Collection Site
code as provided by the Customer - alpha/numeric -->
  <CollectionSiteName>Main Office</CollectionSiteName>
  <CollectionSiteAddress1>123 Main St.</CollectionSiteAddress1>
  <CollectionSiteAddress2>Unit D</CollectionSiteAddress2>
  <CollectionSiteCity>Barrie</CollectionSiteCity>
  <CollectionSiteStateProv>ON</CollectionSiteStateProv>
  <CollectionSiteZipPostalCode>L4N 8W8</CollectionSiteZipPostalCode>
  <CollectionSiteCountry>CA</CollectionSiteCountry>
  <CollectionSitePhone>705-444-2121</CollectionSitePhone>
</CollectionSite>
<Collector>
  <CollectorName>Jennifer Smith</CollectorName>
  <CollectorPhone>705-555-1212</CollectorPhone>
  <CollectorFax>705-555-2121</CollectorFax>
  <CollectorOther>alpha/numeric</CollectorOther> <!-- Text field as entered on CCF -->
  <CollectorEmail>collector@collection.com</CollectorEmail> <!-- if provided by lab from
eCCF -->
</Collector>
<Lab>

```

```

<LabName>Clinical Reference Laboratory, Inc.</LabName>
<LabId>45</LabId>
<LabAddress1>123 Main St.</LabAddress1>
<LabAddress2>Unit D</LabAddress2>
<LabCity>Barrie</LabCity>
<LabStateProv>ON</LabStateProv>
<LabZipPostalCode>L4N 8W8</LabZipPostalCode>
<LabCountry>CA</LabCountry>
<LabPhone>705-444-2121</LabPhone>
<LabFax>705-444-1212</LabFax>
<LabAccountNumber>45678</LabAccountNumber>
<LabSubaccountNumber>9876</LabSubaccountNumber>
<LabSpecimenId>CF18142152</LabSpecimenId> <!-- IMPORTANT SpecimenId is not a unique
value and should not be relied upon as a record reference -->
<LabAccessionNumber>12345678</LabAccessionNumber>
<LabCustomerReferenceId>AU123456</LabCustomerReferenceId> <!-- reference number provided
to the lab, by the customer, on a per test basis-->
<LabCertifyingScientist>Brianna Dingman</LabCertifyingScientist>
</Lab>
<Results>
  <MroVerification>
    <MroName>Dr. Alfred Smith</MroName>
    <MroQualifyingOrganization>MROCC</MroQualifyingOrganization>
    <MroRegistrationNumber>1234567890</MroRegistrationNumber>
    <MroOverallResult>Positive Dilute</MroOverallResult> <!-- MRO account defined result
-->
    <MroMisOverallResult>Positive</MroMisOverallResult> <!-- Mapped DOT MIS result value;
Positive, Negative, Refusal - Substituted, Refusal - Adulterated, Refusal - Shy Bladder,
Other Refusals, Cancelled -->
  </MroVerification>
  <LabPanel> <!-- File may contain multiple lab panels depending on tests ordered and
additional tests performed (e.g. confirmation tests)-->
    <Panel>
      <LabPanelCode>W215</LabPanelCode>
      <LabPanelName>Urine Drug Test</LabPanelName>
      <LabSpecimenType>Urine</LabSpecimenType>
    </Panel>
    <Result>
      <LabTestName>6-Monoacetylmorphine</LabTestName>
      <LabTestCode>6AM</LabTestCode>
      <UnitOfMeasure>ng/mL</UnitOfMeasure>
      <ScreenCutoff>10</ScreenCutoff>
      <ConfirmCutoff>10</ConfirmCutoff>
      <MroResult>Negative</MroResult>
      <NonRegQuant/> <!-- if "Regulated=false" & "EmployerNonRegQuant=true" then a
quantitative, numeric value for the confirmed drug will be included -->
    </Result>
    <Result>
      <TestCodeName>Amphetamines</TestCodeName>
      <TestCodeNumber>AMP</TestCodeNumber>
      <UnitOfMeasure>ng/mL</UnitOfMeasure>
      <ScreenCutoff>500</ScreenCutoff>

```

```

    <ConfirmCutoff>250</ConfirmCutoff>
    <MroResult>Negative</MroResult>
    <NonRegQuant/>
  </Result>
  <Result>
    <TestCodeName>Cocaine</TestCodeName>
    <TestCodeNumber>COC</TestCodeNumber>
    <UnitOfMeasure>ng/mL</UnitOfMeasure>
    <ScreenCutoff>150</ScreenCutoff>
    <ConfirmCutoff>100</ConfirmCutoff>
    <MroResult>Positive</MroResult>
    <NonRegQuant/>
  </Result>
  <Result>
    <TestCodeName>Cannabis</TestCodeName>
    <TestCodeNumber>THC</TestCodeNumber>
    <UnitOfMeasure>ng/mL</UnitOfMeasure>
    <ScreenCutoff>50</ScreenCutoff>
    <ConfirmCutoff>15</ConfirmCutoff>
    <MroResult>Negative</MroResult>
    <NonRegQuant/>
  </Result>
  <Result>
    <TestCodeName>Opiates</TestCodeName>
    <TestCodeNumber>OPI</TestCodeNumber>
    <UnitOfMeasure>ng/mL</UnitOfMeasure>
    <ScreenCutoff>2000</ScreenCutoff>
    <ConfirmCutoff>2000</ConfirmCutoff>
    <MroResult>Negative</MroResult>
    <NonRegQuant/>
  </Result>
  <Result>
    <TestCodeName>Oxycodone</TestCodeName>
    <TestCodeNumber>OXY</TestCodeNumber>
    <UnitOfMeasure>ng/mL</UnitOfMeasure>
    <ScreenCutoff>100</ScreenCutoff>
    <ConfirmCutoff>100</ConfirmCutoff>
    <MroResult>Negative</MroResult>
    <NonRegQuant/>
  </Result>
  <Result>
    <TestCodeName>Phencyclidine</TestCodeName>
    <TestCodeNumber>PCP</TestCodeNumber>
    <UnitOfMeasure>ng/mL</UnitOfMeasure>
    <ScreenCutoff>25</ScreenCutoff>
    <ConfirmCutoff>25</ConfirmCutoff>
    <MroResult>Negative</MroResult>
    <NonRegQuant/>
  </Result>
  <Result>
    <TestCodeName>MDMA/MDA</TestCodeName>
    <TestCodeNumber>MDMA</TestCodeNumber>

```

```

    <UnitOfMeasure>ng/mL</UnitOfMeasure>
    <ScreenCutoff>500</ScreenCutoff>
    <ConfirmCutoff>250</ConfirmCutoff>
    <MroResult>Negative</MroResult>
    <NonRegQuant/>
  </Result>
  <Result>
    <TestCodeName>Hydrocodone</TestCodeName>
    <TestCodeNumber>HCO</TestCodeNumber>
    <UnitOfMeasure>ng/mL</UnitOfMeasure>
    <ScreenCutoff>300</ScreenCutoff>
    <ConfirmCutoff>100</ConfirmCutoff>
    <MroResult>Negative</MroResult>
    <NonRegQuant/>
  </Result>
  <Result>
    <TestCodeName>Creatinine</TestCodeName>
    <TestCodeNumber>CR</TestCodeNumber>
    <UnitOfMeasure>mg/dL</UnitOfMeasure>
    <ScreenCutoff>2-20</ScreenCutoff>
    <ConfirmCutoff/>
    <MroResult>Dilute</MroResult>
  </Result>
  <Result>
    <TestCodeName>Specific Gravity</TestCodeName>
    <TestCodeNumber>SG</TestCodeNumber>
    <UnitOfMeasure/>
    <ScreenCutoff>1.0020-1.0030</ScreenCutoff>
    <ConfirmCutoff/>
    <MroResult>Dilute</MroResult>
  </Result>
</LabPanel>
<LabPanel>
  <Panel>
    <LabPanelCode>CONFIRM</LabPanelCode> <!-- Lab confirmation test code -->
    <LabPanelName>Urine Confirmation Test</LabPanelName>
  </Panel>
  <Result>
    <TestCodeName>Hydrocodone</TestCodeName>
    <TestCodeNumber>HCO</TestCodeNumber>
    <UnitOfMeasure>ng/mL</UnitOfMeasure>
    <ScreenCutoff>300</ScreenCutoff>
    <ConfirmCutoff>100</ConfirmCutoff>
    <MroResult>Positive</MroResult>
    <NonRegQuant/>
  </Result>
</LabPanel>
</Results>
<ReportTimeline>
  <SpecimenCollectedDate>05-23-2024</SpecimenCollectedDate> <!-- Dates in mm-dd-yyyy
format -->
  <MroReceivedCcf2Date>05-23-2024</MroReceivedCcf2Date>

```

```

<LabReceivedDate>05-24-2024</LabReceivedDate>
<LabReceivedTime>17:05:00Z</LabReceivedTime> <!-- Times in HH:mm:ssZ (Zulu/GMT) format
-->
<LabResultedToMroDate>05-26-2024</LabResultedToMroDate>
<LabResultedToMroTime>10:09:05Z</LabResultedToMroTime>
<MroReportedDate>05-26-2024</MroReportedDate>
<MroReportedTime>16:30:09Z</MroReportedTime>
<MroRecalledDate/> <!-- If MRO recalls the report, the date and time (realted to the
"ReportRevisionNumber") recalled date will be noted here. The above "MroReportedDate" will
reflect the revised report's date of reporting -->
<MroRecalledTime/>
<MroRereleasedDate/> <!-- if the MRO report has been recalled, this is the date and time
it was re-released -->
<MroRereleasedTime/>
</ReportTimeline>
<MroNotes>
  <MroLetterComment>MRO letter comment will appear here.</MroLetterComment>
  <MroSafetyAdvisoryComment>MRO Safety Advisory comment will appear
here.</MroSafetyAdvisoryComment>
</MroNotes>
<Documents>
  <MroLetter>Base64 Encoded PDF</MroLetter>
  <SignedCcf2>Base64 Encoded JPG</SignedCcf2>
  <Cg719P>Base64 Encoded PDF</Cg719P>
</Documents>
</VeriportMroData>

```